

Grainger Ultrasonic Cleaning Bath Accessories



General Description

Grainger accessories are designed to maximize the effectiveness of your baths. It is important that parts do not rest on the bottom of the bath, as it can inhibit the ultrasonic activity and possibly damage the unit.

A variety of perforated trays, mesh baskets, beakers, and support racks give you the option of direct or indirect cleaning. All five models come with covers to allow the cleaning bath to heat up more quickly.

Direct Cleaning Method

- Suspends parts using a perforated tray, mesh basket, or a support rack.
- Permits sound waves to penetrate directly to the parts.

Indirect Cleaning Method

- Suspends parts using a solid tray, or beakers with a beaker cover.
- Ideal for cleaning small parts.
- Discarded soil cleans easily from a beaker or solid tray.
- Indirect method is ideal for sample preparation.

Additional Features

- Perforated trays, solid insert trays, and wire mesh baskets have built-in handles specially designed to prevent accessories from touching the tank bottom.
- Support racks allow total volume of the tank to be used while keeping parts, fixtures, or containers about ½" above the tank bottom.
- Support racks also are useful as countertop drain racks, holding parts above the countertop surface.

AVAILABLE THROUGH

GRAINGER

FOR THE ONES WHO GET IT DONE

Grainger – Order using part numbers listed under model numbers

	Ordering Part No.	Model 1800	Model 2800	Model 3800	Model 5800	Model 8800
	Cover (Supplied with Unit)	3KWH1	3KWH7	3KWK4	3KWN4	3KWP4
	Perforated Tray	3KVZ4	3KWC1	3KWD1	3KWD7	3KWE7
	Mesh Basket	3KVZ7	3KWC4	3KWL4	—	—
	Solid Insert Tray	20KN62	3KWJ1	3KWK7	20KN63	3KWP7
	Beaker Positioning Cover	(1x600 ml) 3KWH4	(2x600 ml) 3KWJ4 (2x250 ml) 20KN65	(2x600 ml) 20KN64 (3x250 ml) 3KWL1	(4x600 ml) 3KWN7	(6x600 ml) (SS) 3KWR1
	Support Racks	20KN66	20KN67	20KN68	20KN69	20KN70
	Beakers (All Models)	250 ml Glass – 20KN71 400 ml PP – 20KN72 600 ml Glass – 20KN73 600 ml SS – 20KN74				