



Filter Holders

Glass Microanalysis, 47mm

SIBATA 6168-Series Filter Holders are offered with a wide variety of components for different applications. Available with fritted glass support base for general and biological applications, or with removable stainless steel screen support base and PTFE gasket for ultraclean filtrate or proteinaceous solutions. The 300mL borosilicate glass funnel is graduated from 75 to 300mL in 25mL subdivisions.



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veege.com • email: sales@veege.com

BRC06092014
Printed in the USA



6168-4711
Funnel, 300ml



6168-4751
No 8 Silicon Stopper

6168-4762
Aluminum Clamp
Anodized 47mm I.D.



6168-4741
PTFE Gasket



6168-4721
Glass Support Base
For Stainless Steel Screen

6168-4731
Stainless Steel Screen

Specifications

- Filter Diameter: 47mm
- Pre-filter Diameter: 35mm
- Effective Filtration Area: 9.6cm²
- Maximum Capacity: 350mL
- Graduated Volume: 300mL
- Connection: Silicone Stopper (No. 8)
- Sterilization: Autoclavable (at 121°C / 250°F for 20 minutes)

Cat. No.	Description
6168-4711	Funnel, 47mm, 300mL
6168-4721	Glass Support Base, for Stainless Steel Screen, 47mm
6168-4722	Glass Support Base, Fritted, 160-250µm Pore Size, 47mm
6168-4731	Stainless Steel Screen (SUS316)
6168-4741	PTFE Gasket (pack of 10)
6168-4751	No. 8 Silicone Stopper (fits 1000mL Flask)
6168-4762	Anodized Aluminum Clamp

Package One (Cat. No. 6168-0001)

General & Biological Applications

Funnel (cat. no. 6168-4711)
Glass Support Base (Fritted) (cat. no. 6168-4722)
Silicone Stopper (cat. no. 6168-4751)
Anodized Aluminum Clamp (cat. no. 6168-4762)

Package Two (Cat. No. 6168-0002)

Ultraclean Filtrate or Proteinaceous Solutions

Funnel (cat. no. 6168-4711)
Glass Support Base (SS Screen) (cat. no. 6168-4721)
Stainless Steel Screen (cat. no. 6168-4731)
PTFE Gasket (pk/10) (cat. no. 6168-4741)
Silicone Stopper (cat. no. 6168-4751)
Anodized Aluminum Clamp (cat. no. 6168-4762)

Fritted Glass Support Base



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veege.com • email: sales@veege.com